

OxyCide™ Daily Disinfectant Cleaner

A Non-bleach Disinfectant Concentrate for Daily Cleaning with Efficacy Against C. auris and C. difficile spores







Improve Environmental Outcomes

OxyCide[™] Daily Disinfectant Cleaner delivers sporicidal efficacy but has been designed for daily use. It is an EPA-registered, dilutable concentrate that is effective against *Clostridioides difficile* endospores, *Candida auris* and a broad spectrum of other organisms in a non-bleach formula.

- Proactive sporicide use provides daily defense against *Clostridioides difficile* and *Candida auris*.
- EPA-registered, one-step disinfectant cleaner, virucide and deodorizer with sporicidal activity that delivers complete kill in 5 minutes or less.
- The diluted product, as applied according to label instructions, requires no Personal Protective Equipment (PPE).

Standardize and Simplify

With broad spectrum efficacy and superior material compatibility, OxyCide truly enables product standardization and efficiency improvements. Only one product is needed for daily cleaning, discharges and isolation.

- Enables product standardization and simplification, only one product needed for daily cleaning, discharges and isolation.
- Reduces complexity of training on multiple products.
- Delivers improved material compatibility and dries without leaving any visible residue on surfaces.
- Helps reduce risk and cost of replacing damaged goods. Does not corrode surfaces, damage mattresses or bleach soft goods.
- No issues with absorption or binding like quat-based chemistries.

OxyCide kills the following organisms in 5 minutes or less when used as directed:*

BACTERIAL ENDOSPORES

Clostridioides difficile

BACTERIA

Acinetobacter baumannii Bordetella pertussis Carbapenemase producer Klebsiella pneumonia -KPC(a CRE bacteria) Escherichia coli Extended-Spectrum Beta lactamase producing Escherichia coli (ESBL), Klebsiella pneumoniae MDR Acinetobacter baumannii Methicillin Resistant Staphylococcus aureus (MRSA) Methicillin Resistant Staphylococcus epidermidis (MRSE) Proteus mirabilis Pseudomonas aeruginosa Salmonella enterica Staphylococcus aureus (USA300) Staphylococcus aureus (USA400) Staphylococcus aureus (VISA) Staphylococcus aureus Streptococcus pneumoniae Streptococcus pyogenes Vancomycin Resistant Enterococcus faecalis (VRE)

YEAST

Candida albicans Candida auris

VIRUSES

Human Immunodeficiency Virus Type 1 Hepatitis B Hepatitis C Herpes Simplex Type I Human Coronavirus Influenza A virus Respiratory Syncytial Virus Vaccinia Virus Norovirus Rhinovirus Rotavirus

Superior Material Compatibility

OxyCide Daily Disinfectant Cleaner's active ingredients, hydrogen peroxide and peracetic acid, provide favorable material compatibility that helps minimize surface damage and does not leave residual films or salts behind. Material compatibility tests were conducted in a laboratory setting using controlled application techniques. The images below are of the coupons of various materials, comparing OxyCide to other products used in healthcare settings. (See back page for additional details.)

	Chrome Day 14	Chrome Day 28	Stainless Steel Day 14	Stainless Steel Day 28	Day Zero	
OxyCide Daily Disinfectant Cleaner						-
3M Neutral Quat Disinfectant Cleaner Concentrate		0		in the	Chrome Day O Water	Stainless Steel Day O
Clorox Commercial Solutions Ultra Clorox Germicidal Bleach					Chrome Day 14	Chrome Day 28
Clorox Commercial Solutions Clorox Germicidal Wipes			0		Stainless Steel Day 14	Stainless Steel Day 28
Dispatch Hospital Cleaner with Bleach					Surgical Scrub	s
Dispatch Hospital Cleaner Disinfectant Towels with Bleach					OxyCide (1st Tro	eatment)
Oxivir Five 16 Concentrate		• • •		3	Dispatch Hospital Cleaner Disinfectant Towels w/Bleach (1st Treatment)	
Virex II 256		C. M.			Clorox Bleach 1:10 Dilution (1st Treatment)	

TEST OBJECTIVE:

Evaluate the impact of OxyCide Daily Disinfectant Cleaner and several other disinfectant products on materials commonly found in healthcare settings.

MATERIALS TESTED	PRODUCTS USED
304 Stainless Steel 2"x 2" coupons Plated Chrome 2" x 2" coupons Blue Scrub Shirt (65% polyester/35% cotton blend)	OxyCide Daily Disinfectant Cleaner (Ecolab, EPA reg 1677-237) Neutral Quat Disinfectant Cleaner (3M, EPA reg 473-129-10350) Virex II 256 (Diversey, EPA reg 70627-24) Oxivir Five 16 Concentrate (Diversey, EPA reg 70627-58) Clorox Ultra Germicidal Bleach (Clorox Professional Products EPA reg 67619-8) Dispatch Hospital Cleaner Disinfectant with Bleach (Caltech Industries EPA reg 56392-7) Dispatch Wipe (Caltech Industries EPA reg 56392-8)

METHODS AND RESULTS:

HARD SURFACE TESTING: Test coupons had the above products applied to them via the supplied wipe or product was diluted per label instructions and applied with a saturated microfiber cloth. After each application approximately 1 gram of liquid product remained on each coupon. Coupons were allowed to air dry. After coupons were fully dry, the application cycle was repeated. Coupons were visually assessed for damage and photographs were taken of each coupon. As shown in the pictures on the previous page, the bleach- based products resulted in significantly more visible and irreversible corrosion/discoloring of the coupon materials than did OxyCide Daily Disinfectant Cleaner. Additionally, OxyCide Daily Disinfectant Cleaner did not leave behind a visible residue like the Hydrogen Peroxide or Quat-based products.

SOFT SURFACE TESTING: Approximately 2 grams each of OxyCide Daily Disinfectant Cleaner and a 1:10 dilution of concentrated bleach were applied to separate areas of a blue uniform scrub shirt. A disposable bleach wipe was placed on and pressed against a third area of the scrub shirt. Each spot was allowed to dry and received only one product application. After drying, the shirt was visually assessed for damage and photographs of each area were taken. As shown in the pictures on the previous page, the bleach-based products show irreversible bleaching of the uniform shirt. OxyCide Daily Disinfectant Cleaner showed no visible discoloration of the uniform shirt.



ORDERING INFORMATION: 6000189 2x1 gal (2x96 US fl oz)

WORLD HEADQUARTERS 1 Ecolab Place St. Paul, MN 55102 866 781 8787 www.ecolab.com/oxycide

