

IoShield® (Io-Shield® D)

Double-performing shield to disinfect and protect teats after milking.



- ▲ High disinfectant efficacy of iodine
- ▲ Long lasting Shield protection
- ▲ Easy to apply and remove



APPLICATION

IoShield® is a ready-to-use product. After each milking apply **IoShield®** on the whole teat and do not wipe it. Keep the animals standing or walking for 5 min to ensure biocidal effect.

FEATURES AND BENEFITS

IoShield® contains iodine, which due to its unique PVP carrier eliminates bacteria and disinfects the teats.

The quickly drying, polymerizing Shield of **IoShield®** protects the teats between milking. It is used as part of a herd health plan for the prevention of mammary infections and provides protection from pathogens outside the udder and the environment.

IoShield® provides prevention and protection irrespective of any cow, herd or environmental challenges. Well visible on the teats once applied, **IoShield®** allows a quick and safe milking routine. Easy to remove, the formulation enables easier cleaning.



SPECIFICATIONS

IoShield®	
Physical form	Dark brown liquid
Active substance	Iodine (PVP)
pH at 20°C	5.0 - 5.5
Density at 20°C	0.99 - 1.03
Storage temperature	+5°C to +25°C

LOGISTICS

Packaging	10 kg	20 kg	58 kg	210 kg	1000 kg
Product code	2288660	2267840	2267880	2267920	2288650
EAN code	4028162288667	4028162267846	4028162267884	4028162267921	4028162288650

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtstrasse 7
CH-8304 Wallisellen
www.ecolab.com